

StateMix

Mix With Authority!



Mix and degass hard to blend materials with no blades over a wide range of viscosities in seconds!

- Polyurethanes
- Silicones
- Inks
- Waxes
- Filled polymers
- Epoxies
- And much, much more!



StateMix

1390 Church Ave Winnipeg, Manitoba R2X 1G4, Canada

Ph: (204) 697-2586 Fax: (204) 697-2592

Email: state@mts.net

www.statemix.com

Vortex Mixer

VM-20,000
COMING
SOON!

VM-200

A bench-top lab machine that mixes and degasses up to 200g of material in a plastic cup. This unit is ideal for material labs and quality control applications.

- 3/4 hp, Inverter Duty, Single Phase, 50-60 Cycle, 110-120 Volt or 208-240 Volt
- Mixing Capacity 30 grams to 200 grams
- Adjustable counterweight system
- Active dampening technology
- Frequency Inverter
- Vibration Sensor
- Dynamic Braking System
- CE Certified

DIMENSIONS: 21" TALL 16.3" WIDE 16.3" DEEP



VM-1000

A stand-alone mixer that mixes and degasses up to 1000 grams or 2.2 lbs of material per batch. This unit is used in both lab and for small production runs.

- 2 hp, Inverter Duty, Single Phase, 50-60 Cycle, 208-240 Volt
- Mixing Capacity 300 grams to 1000 grams
- Adjustable counterweight system
- Active dampening technology
- Frequency Inverter
- Vibration Sensor
- Dynamic Braking System
- CE Certified

DIMENSIONS: 37" TALL 27" WIDE 27" DEEP



VM-2500

Developed for larger batch production runs using a variety of materials. The VM-2500 is capable of mixing and degassing up to 2.5 kg or 5.5 lbs of material in a standard 1 gallon or 4 L pail.

- 5 hp, Inverter Duty, 3 Phase, 50-60 Cycle, 208-240 Volt or 400-480 Volt
- Mixing Capacity 500 grams to 2500 grams
- Adjustable counterweight system
- Active dampening technology
- Frequency Inverter
- Vibration Sensor
- Dynamic Braking System
- CE Certified

DIMENSIONS: 40" TALL 36" WIDE 36" DEEP



VM-10,000

Developed for large production runs and is capable of mixing and degassing up to 10 kg or 22 lbs of material per batch in a standard 5 gallon or 20 L plastic pail.

- 10 hp, Inverter Duty, 3 Phase, 50-60 Cycle, 208-240 Volt or 400-480 Volt
- Mixing Capacity 3 kg to 10 kg
- Adjustable counterweight system
- Active dampening technology
- Frequency Inverter
- Vibration Sensor
- Dynamic Braking System
- CE Certified

DIMENSIONS: 48" TALL 48" WIDE 48" DEEP



STATE MIX BUSINESS PARTNER



MITEX GmbH · Ludenberger Str. 36-40 · D-40699 Erkrath
Tel.: +49 (0)211 90006-0 · Fax: +49 (0)211 246118
info@mitex.de · www.mitex.de · www.vortex-mixer.de